
EN Medical device **CZ** Zdravotnické prostředky
DE Medizinprodukt **DK** Medicinsk udstyr **ES**
Dispositivo médico **FI** Lääketieteellinen laite **FR**
Dispositif Medical **GR** Ιατρικοτεχνολογική συσκευή
IS Lækningsbúnaður **IT** Dispositivo medico **JP** 医療
機器 **NL** Medische toepassing **NO** Medisinsk enhet
SE Medicinsk utrustning **SI** Medicinski pripomoček
AE جهاز طبي

EN Do not re-use **CZ** Nepoužívejte znovu **DE**
Nicht wiederverwenden **DK** Må ikke genbruges **ES**
No reutilizar **FI** Ei saa käyttää uudelleen **FR** Usage
unique **GR** Μην επαναχρησιμοποιείτε **IS** Endurnotið
ekki **IT** Non riutilizzabile **JP** 再使用禁止 **NL** Niet
hergebruiken **NO** Ikke for gjenbruk **SE** Återanvänd
ej **SI** Ni namenjeno ponovni uporabi **AE** لا يعاد استخدامه

EN Temperature limit **CZ** Teplotní limit **DE**
Lagerungstemperatur **DK** Temperaturgrænse
ES Temperatura límite **FI** Lämpötilaraja **FR**
Temperature de conservation **GR** Όριο
θερμοκρασίας **IS** Hitamörk **IT** Limite superiore di
temperatura **JP** 温度制限 **NL** Temperatuurgrens
NO Temperaturgrense **SE** Max temperatur **SI**
Temperaturna omejitev **AE** حد درجة الحرارة

EN Keep dry **CZ** Udržujte v suchu **DE** Trocken
lagern **DK** Opbevares tørt **ES** Mantener seco
FI Pidä kuivana **FR** Conserver à l'abri de
l'humidité **GR** Διατηρείτε στεγνό **IS** Haldið þurru **IT**
Mantenere asciutto **JP** 湿気厳禁 **NL** Droog houden
NO Holdes tørr **SE** Förvaras torr **SI** Hranite na
suhem mestu **AE** يحفظ جافاً

EN Capable of being recycled (carton only)
CZ lze recyklovat (pouze kartony) **DE** Karton
kann recycelt werden **DK** Kan genanvendes
(kun kartonen) **ES** Se puede reciclar **FI**
Voidaan kierrättää (vain pahvi) **FR** emballage
carton recyclable **GR** Αβακυκλώσιμο (μόνο
χαρτοκιβώτιο) **IS** Endurvinnanlegt (aðeins öskju)
IT Riciclabile (solo cartone) **JP** 外箱のみリサイクル
可能 **NL** Geschikt voor recycling (verpakking) **NO**
Kan resirkuleres (bare kartong) **SE** Kan återvinnas
(endast kartong) **SI** Primerno za reciklažo (samo
škafka) **AE** قابل لإعادة التدوير

eakin® pouches/ eakin® ‘pelican pouches™’: flat and soft convex ostomy pouches

ABOUT THIS PRODUCT.

eakin® pouches are used to collect output from a stoma. The baseplate attached to the pouch is intended to adhere to peristomal skin.

Individual wear times may vary.

eakin® soft convex pouches may be suitable for retracted, flush or other challenging stomas only after assessment from a stoma care professional.

eakin® pre-cut pouches may be suitable for stomas with a regular shape and size. **eakin®** cut-to-fit pouches allow a customised fit.

For those requiring additional security a stoma belt may be used. (Product availability may vary in your country).

WARNINGS AND CAUTIONS

Consult a healthcare professional if you have any skin irritation, redness or discomfort in the area of product application prior to use. If during use you experience any skin irritation, redness or discomfort in the area of product application, discontinue use and consult a healthcare professional. Do not apply to broken skin or open wounds.

For single use only. Re-use of a single use product may have potential to cause harm to the user due to the risk of contamination and impairment in product function.

Skin adhesion and wear time may be reduced in conditions of elevated temperature or humidity, or if the baseplate is re-positioned after first application.

OTHER INFORMATION

Store between 10 and 30 degrees Celsius. Store flat in the carton. Do not refrigerate.

eakin® accepts no liability for any injury or other loss that may arise if this product is used in a manner contrary to **eakin®**'s current recommendations.

Any serious incident encountered, in relation to the use of this device, should be reported to the manufacturer and the competent authority in the country of use. A serious incident is one involving serious deterioration in health or patient death.

Instructions for use

1. PREPARATION

Ensure skin is clean and dry before application (1.a).

Use the enclosed stoma size measuring guide to measure the size and shape of the stoma (1.b). It is important to ensure an exact fit between the baseplate hole and the stoma to reduce the risk of leaks and skin problems.

If your stoma is not an exact circle or oval as shown on the stoma size measuring guide, you will need to follow the instructions on the enclosed stoma size cutter guide to size your stoma and your pouch.

Using the cutting guide on the baseplate, cut a hole to fit the exact pre-measured size and shape of the stoma using curved scissors, taking care not to cut through the pouch film on the other side (1.c).

For pre-cut pouches select the variant which matches the size of the stoma.

2. APPLICATION

Peel off the protective film from the baseplate (2.a).

Align the centre hole of the baseplate around the stoma.

Starting at the bottom of the baseplate press gently against your skin, smoothing upwards and outwards to ensure the entire baseplate is securely adhered to the skin. (2.b).

(Filter covers are included which can be used when bathing or to prevent 'pancaking'. Apply the cover by lifting up the split fabric on the pouch and placing the cover over the filter which can be seen as a dark circle).

3.1 EMPTYING DRAINABLE POUCHES

Peel both hook tabs off the velcro strip (3.a). Unfold the outlet in the position where you want to empty the pouch. Hold the bottom of the outlet with one hand and pull upwards on the opening tab to make a larger opening for emptying (3.b). Allow contents to drain into the toilet. Clean the outlet in one sweeping movement away from your body.

To close: fold the outlet back up and secure by attaching the two hook tabs onto the velcro strip. Tuck the velcro strip up under the split at the bottom of the pouch fabric and attach to the velcro dot for added comfort and discretion (3.c). (Tuck-up option only available on bags with fabric on 'non-body' side).

3.2 EMPTYING UROSTOMY POUCHES

Pull the bung closure out of the bung to allow emptying.

To control the flow of urine from the bung, simply squeeze. To close, re-insert the bung closure.

4. REMOVAL

Gently remove the pouch by peeling the baseplate downwards and away from the skin. Apply light pressure to your skin with your free hand to aid removal (4.a).

5. DISPOSAL

Dispose of waste product in accordance with local guidelines. Do not flush the product down the toilet.

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© 2020 TG Eakin Ltd.
TG Eakin Limited,