MD **EN** Medical device **CZ** Zdravotnické prostředky
DE Medizinprodukt **DK** Medicinsk udstyr
ES Dispositivo médico **FI** Lääketieteellinen laite
FR Dispositif Médical **GR** Ιατρικότεχνολογική
 συσκευή **IS** Lækningsgáðunáður **IT** Dispositivo
 medico **JP** 医療機器 **NL** Medisch hulpmiddel
NO Medisinsk enhet **SE** Medicinsk utrustning
SI Medicinski pripomoček **AE** جهاز طبي

 **EN** Do not re-use **CZ** Nepoužívejte znovu **DE** Nicht
 wiederverwenden **DK** Må ikke genbruges **ES** No
 reutilizar **FI** Ei saa käyttää uudelleen **FR** Usage
 unique **GR** Μην επαναχρησιμοποιείτε **IS** Endurnotit
 ekki **IT** Non riutilizzabile **JP** 再使用禁止 **NL** Niet
 hergebruiken **NO** Ikke for gjenbruk **SE** Återanvänd
 ej **SI** Ni namenjeno ponovni uporabi **AE** لا يُعاد استخدامه

 **EN** Temperature limit **CZ** Teplotní limit
DE Lagerungstemperatur **DK** Temperaturgrænse
ES Temperatura límite **FI** Lämpötilaraja
FR Température de conservation
GR Όριο θερμοκρασίας **IS** Hitamörk
IT Limite superiore di temperatura **JP** 温度制限
NL Bewaarttemperatuur **NO** Temperaturgrense
SE Max temperatur **SI** Temperaturna omejitev
AE حد درجة الحرارة

 **EN** Keep dry **CZ** Udržujte v suchu **DE** Trocken
 lagern **DK** Opbevares tørt **ES** Mantener seco
FI Pidä kuivana **FR** Conserver à l'abri de
 l'humidité **GR** Διατηρείται στεγνό **IS** Haldið þurru
IT Mantenere asciutto **JP** 湿気厳禁 **NL** Droog
 bewaren **NO** Holdes tørt **SE** Förvaras torrt
SI Hranite na suhem mestu **AE** يحفظ جافاً

 **EN** Capable of being recycled (carton only)
CZ lze recyklovat (pouze kartony)
DE Karton kann recycelt werden
DK Kan genanvendes (kun kartoner)
ES Se puede reciclar
FI Voidaan kierrättää (vain pahvi)
FR emballage carton recyclable
GR Ανακυκλώσιμο (μόνο χαρτοκιβώτιο)
IS Endurninnanlegt (aðeins óskju)
IT Riciclabile (solo cartone)
JP 外箱のみリサイクル可能
NL Geschikt voor recycling (verpakking)
NO Kan resirkuleres (bare kartong)
SE Kan återvinnas (endast kartong)
SI Primerno za reciklažo (samo škatla)
AE قابل لإعادة التدوير

CH IMPORT

Publicare AG
 Vorder Böde 9
 3452 Oberrohrdorf
 Switzerland

EN For Switzerland only **CZ** Pouze pro použití
 ve Švýcarsku **DE** Nur gültig in der Schweiz
DK Kun til Schweiz **ES** Solo aplicable a residentes de
 Suiza **FI** Koskee ainoastaan Sveitsiä **FR** Pour utilisation en
 Suisse uniquement **GR** Μόνο για χρήση στην Ελλάδα
IS Eingöngu fyrir Sviss **IT** Solo per la Svizzera
JP スイス居住のユーザー対象 **NL** Alleen voor gebruik in
 Zwitserland **NO** Bare for brukere i Sveits **SE** Gäller endast
 Schweiz **SI** Samo za uporabo v Švici
AE موقوف لنيروسيويسلأى

eakin® pouches/ eakin® ‘pelican pouches™’: flat and soft convex ostomy pouches

ABOUT THIS PRODUCT.

eakin® pouches are used to collect output from a stoma. The baseplate attached to the pouch is intended to adhere to peristomal skin.

Individual wear times may vary.

eakin® soft convex pouches may be suitable for retracted, flush or other challenging stomas only after assessment from a stoma care professional.

eakin® pre-cut pouches may be suitable for stomas with a regular shape and size. **eakin®** cut-to-fit pouches allow a customised fit.

For those requiring additional security a stoma belt may be used. (Product availability may vary in your country).

Intended users - For ostomates, healthcare professionals and caregivers.

WARNINGS AND CAUTIONS

Consult a healthcare professional if you have any skin irritation, redness or discomfort in the area of product application prior to use. If during use you experience any skin irritation, redness or discomfort in the area of product application, discontinue use and consult a healthcare professional. Do not apply to broken skin or open wounds.

For single use only. Re-use of a single use product may have potential to cause harm to the user due to the risk of contamination and impairment in product function.

Skin adhesion and wear time may be reduced in conditions of elevated temperature or humidity, or if the baseplate is re-positioned after first application.

OTHER INFORMATION

Store between 10 and 30 degrees Celsius. Store flat in the carton. Do not refrigerate.

eakin® accepts no liability for any injury or other loss that may arise if this product is used in a manner contrary to **eakin®**'s current recommendations.

Any serious incident encountered, in relation to the use of this device, should be reported to the manufacturer and the competent authority in the country of use. A serious incident is one involving serious deterioration in health or patient death.

Instructions for use

1. PREPARATION

Ensure skin is clean and dry before application (1.a).

Use the enclosed stoma size measuring guide to measure the size and shape of the stoma (1.b). It is important to ensure an exact fit between the baseplate hole and the stoma to reduce the risk of leaks and skin problems.

If your stoma is not an exact circle or oval as shown on the stoma size measuring guide, you will need to follow the instructions on the enclosed stoma size cutter guide to size your stoma and your pouch.

Using the cutting guide on the baseplate, cut a hole to fit the exact pre-measured size and shape of the stoma using curved scissors, taking care not to cut through the pouch film on the other side (1.c).

For pre-cut pouches select the variant which matches the size of the stoma.


2. APPLICATION

Peel off the protective film from the baseplate (2.a).

Align the centre hole of the baseplate around the stoma.

Starting at the bottom of the baseplate press gently against your skin, smoothing upwards and outwards to ensure the entire baseplate is securely adhered to the skin. (2.b).

(Filter covers are included which can be used when bathing or to prevent 'pancaking'. Apply the cover by lifting up the split fabric on the pouch and placing the cover over the filter which can be seen as a dark circle).

3.1 EMPTYING DRAINABLE POUCHES

Peel both hook tabs off the velcro strip (3.a). Unfold the outlet in the position where you want to empty the pouch. Hold the bottom of the outlet with one hand and pull upwards on the opening tab to make a larger opening for emptying (3.b). Allow contents to drain into the toilet. Clean the outlet in one sweeping movement away from your body.

To close: fold the outlet back up and secure by attaching the two hook tabs onto the velcro strip. Tuck the velcro strip up under the split at the bottom of the pouch fabric and attach to the velcro dot for added comfort and discretion (3.c). (Tuck-up option only available on bags with fabric on 'non-body' side).

3.2 EMPTYING UROSTOMY POUCHES

Pull the bung closure out of the bung to allow emptying.


To control the flow of urine from the bung, simply squeeze. To close, re-insert the bung closure.

4. REMOVAL

Gently remove the pouch by peeling the baseplate downwards and away from the skin. Apply light pressure to your skin with your free hand to aid removal (4.a).

5. DISPOSAL

Dispose of waste product in accordance with local guidelines. Do not flush the product down the toilet.

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© 2021 TG Eakin Ltd.
TG Eakin Limited,